

Automated Chemistry Efforts at LLNL Using Eichrom's DGA Resin

Presented to:
SEVENTH WORKSHOP ON THE CHEMISTRY OF
THE HEAVIEST ELEMENTS
October 13, 2009

**Roger A. Henderson, Kenton J. Moody, Dawn A. Shaughnessy,
Megan Bennett*, and Julie M. Gostic***

Chemistry, Materials, Earth, and Life Sciences Directorate
Lawrence Livermore National Laboratory


*University of Nevada-Las Vegas

This work was performed under the auspices of the U.S. Department of Energy by
Lawrence Livermore National Laboratory under Contract DE-AC52-07NA27344.
Lawrence Livermore National Laboratory, P.O. Box 808, Livermore, CA 94551


LLNL-PRES-406833

Automated Chemistry Efforts at LLNL Using Eichrom's DGA Resin

- Chemistry Desires for Automation
 - Speed
 - Reproducibility
 - Lower Dose to Operators
- Begin with a Heavy Element application
 - Single Column chemistry
 - Automated column regeneration
 - Long operating periods (months!)
- Basic Design



Automated Chemistry Efforts at LLNL Using Eichrom's DGA Resin




Automated Chemistry Efforts at LLNL Using Eichrom's DGA Resin


- System Components
 - Cheminert valves, PEEK material, micro-electric actuation
 - ISCO peristaltic pumps
 - Tubing – FEP – ID = 0.03” (0.762 mm)
 - Columns – FEP tubing, 10cm length
 - Volumes
 - 0.762 mm tubing – 45.6 μ l
 - 0.508 mm tubing – 20.3 μ l
 - 0.254 mm tubing – 5.1 μ l
 - Frits
 - PEEK, 0.5 μ or 2 μ Frit-In-A-Ferrule™



Automated Chemistry Efforts at LLNL Using Eichrom's DGA Resin



Automated Chemistry Efforts at LLNL Using Eichrom's DGA Resin




Automated Chemistry Efforts at LLNL Using Eichrom's DGA Resin

- System Notes
 - Resin slurry concentration is key
 - Slurry tends to settle out of the load tube from the reservoir
 - DGA seems to work OK; AG-1 does not behave as well
 - As the column fills system back pressure increases




Automated Chemistry Efforts at LLNL Using Eichrom's DGA Resin



DGA
Cartridge
Column




Automated Chemistry Efforts at LLNL Using Eichrom's DGA Resin



Automated Chemistry Efforts at LLNL Using Eichrom's DGA Resin

Test 3



Automated Chemistry Efforts at LLNL Using Eichrom's DGA Resin

- Items for Study
 - System dead volumes, how they relate to elution parameters
 - Off-line vs. On-line operation
 - Couple now with automated column generation
 - Kinetics – how fast can we run the pumps?
 - Chemical – Affect on separation capability
 - Mechanical – Column support frit size



Automated Chemistry Efforts at LLNL Using Eichrom's DGA Resin

- Future Plans
 - Sample Introduction System
 - Gas jet into a frit
 - Continuous deposition into load solution?
 - Sample Detection System
 - Functionalized surfaces on glass
 - Fast Kinetics
 - Additional Selectivity?
 - On-line tests at an accelerator facility
 - Carrier free activity (not the case currently with Ta)
- Questions?

